

**INTERNATIONAL
PATENT**

Item Code 910004



FlowMIR[®]

Disposable Turbine Flowmeter
MIR exclusive product
International patent

Each turbine, which includes a cardboard mouthpiece, has been individually factory tested with a computerized system

For each patient, after the Spirometry test, both turbine and mouthpiece are thrown away



FlowMIR Dispensers (10 and 60 pcs)



Comfortable
Packaging



Singularly Tested
and Packed



Eliminates staff
clean-up time



No Sterilisation



No Calibration



No Cross
Contamination



No Antibacterial
Filter



Always 100%
accurate and hygienic



ATS/ERS Compliant
for accuracy



Not affected by
Vapour Condensation



Not affected by
Ambient conditions



The Best Sensor
for Spirometry