

six⁰⁰ series of Aneroid Sphygmomanometers

ACCOSON

Established 1859
Accuracy starts in the name

The new **six⁰⁰** series is a premium range of aneroid sphygmomanometers offering the clinician the highest quality at an affordable price. A modern solution for all locations with innovative and practical features.

- Made in England and guaranteed accurate to ISO 81060-1
- Range 0 - 300 mmHg, with the scale calibrations starting at zero
- Clear white face with fine red pointer gives pin point visibility
- Smooth lines for easy cleaning without awkward areas
- Designed to withstand double full scale overpressure
- Five year calibration warranty

Optional Wall mounting
velcro[®] plate
Code: 1660



Wall model
Code: 0642



Desk Model
Code: 0632



- Spacious storage for latex free inflation system
- Desk model with practical handle and scratch free feet
- Wall model with mounts for fixed or swivel operation
- Hole free wall mounting available using optional velcro[®] plate
- Mobile stand model with adjustable height and inclined gauge
- Special mounts for rail and other systems on request



Stand model
Code: 0662